

Rabbit Anti-Human Cd56/Ncam1 Monoclonal Antibody For Flow Cytometry - P13591

Item Number: CM0004450

Introduction

Rabbit monoclonal antibody against human CD56/NCAM1, validated for flow cytometry. Recombinant immunogen targeting amino acids 20-413. Reacts with human. Ideal for studying neural cell adhesion and virus-host interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD56/NCAM1 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
GeneID	4684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-413 of human CD56/NCAM1 (NP_851996.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag