


Rabbit Anti-Human Cd58 (Lfa-3) Monoclonal Antibody For Fc - P19256

Item Number: CM0008484



Introduction

-conjugated rabbit monoclonal antibody targeting human CD58 (LFA-3), validated for flow cytometry. Explore T-cell interactions and immune synapse research with this high-specificity reagent.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD58 mAb
Remarks/Alias	CD58; LFA-3; LFA3; ag3; CD58 molecule
Species	Human
Gene ID	965
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-215 of human CD58 (NP_001770.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P19256
Protein Weight	26kDa/27kDa/28kDa
Shipping	Ice bag