


Rabbit Anti-Human Cd58 Mab For Fc - P19256

Item Number: CM0011454



Introduction

Rabbit anti-Human CD58 monoclonal antibody (LFA-3, ag3) validated for flow cytometry with human samples. Detects 28kDa lymphocyte function-associated antigen 3, a CD2 ligand critical for T-cell interactions. Ideal for immune synapse studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD58 mAb
Remarks/Alias	ag3; LFA3; LFA-3
Species	Human
GeneID (Human)	965
GeneID	965
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-215 of human CD58 (NP_001770.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P19256
Protein Weight	28kDa
Shipping	Ice bag