


Rabbit Anti-Human Cd58 Monoclonal Antibody For Fc - P19256

Item Number: CM0005123



Introduction

-conjugated rabbit monoclonal antibody targeting human CD58 (LFA-3). Validated for flow cytometry, ideal for investigating T-cell adhesion, immune synapse dynamics, and lymphocyte interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD58 mAb
Remarks/Alias	ag3; LFA3; LFA-3
Species	Human
Gene ID (Human)	965
Gene ID	965
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-215 of human CD58 (NP_001770.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	FC
Cross Reactivity	Human
SWISS	P19256
Protein Weight	28kDa
Shipping	Ice bag