


Rabbit Anti-Human Cd59 Monoclonal Antibody For Fc - P13987

Item Number: CM0010662



Introduction

Rabbit anti-Human CD59 monoclonal antibody validated for flow cytometry. Targets CD59, a potent MAC inhibitor protecting cells from complement lysis. Ideal for immunology and complement research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD59 mAb
Remarks/Alias	1F5; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; M1RL; HRF20; MACIF; MEM43; MIC11; MSK21; 16.3A5; HRF-20; MAC-IP; p18-20
Species	Human
GeneID (human)	966
GeneID	966
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 26-101 of human CD59 (NP_000602.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13987
Protein Weight	14kDa
Shipping	Ice bag