


Rabbit Anti-Human Cd59 Monoclonal Antibody For Flow Cytometry - P13987

Item Number: CM0011638



Introduction

-conjugated rabbit monoclonal antibody targeting human CD59 (protectin). Recombinant immunization ensures lot-to-lot consistency. Validated for flow cytometry. Suitable for complement research and in vitro diagnostic development. Human cross-reactivity.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD59 mAb
Remarks/Alias	1F5; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; M1RL; HRF20; MACIF; MEM43; MIC11; MSK21; 16.3A5; HRF-20; MAC-IP; p18-20
Species	Human
GeneID (Human)	966
GeneID	966
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-101 of human CD59 (NP_000602.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13987
Protein Weight	14kDa
Shipping	Ice bag