


Rabbit Anti-Human Cd6 Mab For Fc - P30203

Item Number: CM0010733



Introduction

Rabbit anti-Human CD6 monoclonal antibody specifically detects human CD6/TP120 in flow cytometry. Useful for studying T-cell activation, immunological synapse formation, and inflammatory responses, covering isoforms from 58 to 71 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD6 mAb
Remarks/Alias	CD6; TP120; CD6 molecule
Species	Human
Gene ID (Human)	923
Gene ID	923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-364 of human CD6 (NP_006716.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30203
Protein Weight	58kDa/60kDa/63kDa/64kDa/68kDa/71kDa
Shipping	Ice bag