


Rabbit Anti-Human Cd6 Monoclonal Antibody For Fc - P30203

Item Number: CM0005727



Introduction

Rabbit anti-Human CD6 Monoclonal Antibody for FC. Targets T-cell differentiation antigen CD6, involved in T-cell activation and bacterial pattern recognition. Suitable for flow cytometry applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD6 mAb
Remarks/Alias	CD6; TP120; CD6 molecule
Species	Human
Gene ID (Human)	923
Gene ID	923
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-364 of human CD6 (NP_006716.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30203
Protein Weight	58kDa/60kDa/63kDa/64kDa/68kDa/71kDa
Shipping	Ice bag