

Rabbit Anti-Human Cd6 Monoclonal Antibody For Flow Cytometry (Fc) - P30203

Item Number: CM0011008



Introduction

-conjugated rabbit monoclonal antibody targeting human CD6 (T12/TP120). Suitable for flow cytometry (FC), this antibody recognizes the extracellular domain of CD6, a T-cell adhesion and immune regulatory protein. Cross-reacts with human CD6.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD6 mAb
Remarks/Alias	CD6; TP120; CD6 molecule
Species	Human
GeneID (Human)	923
GeneID	923
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-364 of human CD6 (NP_006716.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30203
Protein Weight	58kDa/60kDa/63kDa/64kDa/68kDa/71kDa
Shipping	Ice bag