

Rabbit Anti-Human Cd61 Monoclonal Antibody For Fc - P05106

Item Number: CM0009246



Introduction

-conjugated rabbit monoclonal antibody targeting human CD61 (Integrin beta-3, GPIIIa), validated for flow cytometry. Ideal for studying platelet function, cell adhesion, and integrin signaling in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD61 mAb
Remarks/Alias	GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24
Species	Human
GeneID (Human)	3690
GeneID	3690
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 21-718 of human CD61 (NP_000203.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05106
Shipping	Ice bag