

Rabbit Anti-Human CD61 Monoclonal Antibody For Flow Cytometry - P05106

Item Number: CM0008331



Introduction

-conjugated rabbit monoclonal antibody targeting human CD61 (Integrin beta-3/GPIIIa) for flow cytometry. Validated for FC, enabling studies of platelet aggregation, integrin signaling, and cell adhesion. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD61 mAb
Remarks/Alias	GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24
Species	Human
GeneID (Human)	3690
GeneID	3690
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-718 of human CD61 (NP_000203.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05106
Shipping	Ice bag