


Rabbit Anti-Human Cd62E/E-Selectin Monoclonal Antibody For Fc - P16581

Item Number: CM0008360



Introduction

-conjugated rabbit monoclonal antibody targeting human E-selectin (CD62E). Validated for flow cytometry, reacts specifically with human. Immunogen: recombinant protein covering amino acids 22-230. Ideal for studies on leukocyte adhesion and endothelial biology.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD62E/E-Selectin mAb
Remarks/Alias	SELE; CD62E; ELAM; ELAM1; ESEL; LECAM2; E-selectin
Species	Human
Gene ID	6401 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of human CD62E/E-Selectin (NP_000441.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16581
Protein Weight	66kDa
Shipping	Ice bag