


Rabbit Anti-Human Cd62E/E-Selectin Monoclonal Antibody For Flow Cytometry - P16581

Item Number: CM0010470



Introduction

High-quality -conjugated rabbit anti-human CD62E/E-Selectin monoclonal antibody for flow cytometry. Recognizes human E-selectin extracellular domain (aa 22-230). Ideal for immunoadhesion and endothelial research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD62E/E-Selectin mAb
Remarks/Alias	SELE; CD62E; ELAM; ELAM1; ESEL; LECAM2; E-selectin
Species	Human
GeneID (Human)	6401
GeneID	6401
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of human CD62E/E-Selectin (NP_000441.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16581
Protein Weight	66kDa
Shipping	Ice bag