

Rabbit Anti-Human Cd62L/L-Selectin Monoclonal Antibody For Flow Cytometry - P14151

Item Number: CM0010801



Introduction

-conjugated rabbit monoclonal antibody targeting human CD62L/L-Selectin for flow cytometry. Validated for human cross-reactivity, suitable for studying lymphocyte homing, leukocyte rolling, and cell adhesion in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD62L/L-Selectin mAb
Remarks/Alias	TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1
Species	Human
Gene ID (Human)	6402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14151
Protein Weight	42kDa/43kDa
Shipping	Ice bag