

## Rabbit Anti-Human Cd62P/P-Selectin Monoclonal Antibody For Flow Cytometry - P16109

Item Number: CM0012489



### Introduction

-conjugated rabbit monoclonal antibody targeting human P-selectin (CD62P) for flow cytometry. Validated for human samples, provides specific detection of activated platelets and endothelial cells.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD62P/P-selectin mAb                                                                                       |
| Remarks/Alias    | CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM                                                                              |
| Species          | Human                                                                                                                        |
| GeneID (Human)   | 6403                                                                                                                         |
| GeneID           | 6403                                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 42-771 of human CD62P/P-selectin (NP_002996.2) |
| Source           | Rabbit                                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                                        |
| Application      | FC                                                                                                                           |
| Cross Reactivity | Human                                                                                                                        |
| SWISS            | P16109                                                                                                                       |
| Protein Weight   | 91kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |