


Rabbit Anti-Human Cd63 Monoclonal Antibody For Flow Cytometry - P08962

Item Number: CM0010545



Introduction

-conjugated rabbit monoclonal antibody targeting human CD63 (SWISS P08962), validated for flow cytometry. Ideal for research on tetraspanin-mediated signaling, lysosomal trafficking, and melanoma.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD63 mAb
Remarks/Alias	MLA1; ME491; LAMP-3; OMA81H; TSPAN30
Species	Human
GeneID (human)	967
GeneID	967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human CD63 (NP_001771.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08962
Protein Weight	26kDa
Shipping	Ice bag