

Rabbit Anti-Human Cd66A/Ceacam1 Monoclonal Antibody For Fc - P13688

Item Number: CM0005055



Introduction

-conjugated rabbit anti-human CEACAM1 monoclonal antibody for flow cytometry. Recognizes human CEACAM1 (SWISS: P13688), an immune checkpoint, adhesion molecule, and CEA family member. Suitable for T cell, NK, and epithelial cell analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD66a/CEACAM1 mAb
Remarks/Alias	CEACAM1; BGP; BGP1; BGPI
Species	Human
GeneID (Human)	634
GeneID	634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-428 of human CD66a/CEACAM1 (NP_001703.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13688
Protein Weight	58kDa
Shipping	Ice bag