


Rabbit Anti-Human Cd66A/Ceacam1 Monoclonal Antibody For Fc - P13688

Item Number: CM0004499



Introduction

-conjugated rabbit monoclonal antibody against human CD66a/CEACAM1. Validated for flow cytometry (FC), reacts with human. Recombinant protein immunogen, target 58 kDa. Ideal for immune regulation, angiogenesis, and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD66a/CEACAM1 mAb
Remarks/Alias	CEACAM1; BGP; BGP1; BGPI; carcinoembryonic antigen-related cell adhesion molecule 1
Species	Human
GeneID (Human)	634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-428 of human CD66a/CEACAM1 (NP_001703.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13688
Protein Weight	58kDa
Shipping	Ice bag