


Rabbit Anti-Human Cd66B Monoclonal Antibody For Flow Cytometry - P31997

Item Number: CM0012177



Introduction

-conjugated rabbit monoclonal antibody against human CD66b (CEACAM8, CD67, CGM6). Suitable for flow cytometry. Used in cell adhesion research, especially neutrophil heterophilic interactions. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD66b mAb
Remarks/Alias	CD67; CGM6; CD66b; NCA-95
Species	Human
GeneID (Human)	1088
GeneID	1088
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-319 of human CD66b (NP_001807.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P31997
Protein Weight	38kDa
Shipping	Ice bag