

Rabbit Anti-Human Cd68 Mab - P34810

Item Number: CM0005409

Introduction

-conjugated rabbit monoclonal anti-human CD68 antibody for intracellular flow cytometry. Specific for human CD68 (GP110, LAMP4). Ideal for macrophage and tumor cell dissemination studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD68 mAb
Remarks/Alias	GP110; LAMP4; SCARD1
Species	Human
GeneID	968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-319 of human CD68 (NP_001242.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P34810
Protein Weight	37kDa
Shipping	Ice bag