

Rabbit Anti-Human Cd68 Monoclonal Antibody For Flow Cytometry (Intracellular) - P34810

Item Number: CM0008656

Introduction

-conjugated rabbit monoclonal antibody targeting human CD68 (macrosialin) for intracellular flow cytometry. Detects a 37 kDa transmembrane glycoprotein on macrophages and tumor cells. Ideal for immune profiling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD68 mAb
Remarks/Alias	GP110; LAMP4; SCARD1
Species	Human
GeneID (Human)	968
GeneID	968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-319 of human CD68 (NP_001242.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P34810
Protein Weight	37kDa
Shipping	Ice bag