

Rabbit Anti-Human Cd7 Mab - P09564

Item Number: CM0009750

Introduction

-conjugated rabbit monoclonal antibody against human CD7, validated for flow cytometry. Ideal for studying T-cell and NK cell costimulation, cytokine production, and immune cell profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD7 mAb
Remarks/Alias	GP40, TP41, Tp40, LEU-9
Species	Human
Gene ID (Human)	924
Gene ID	924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag