

Rabbit Anti-Human Cd7 Monoclonal Antibody For Fc - P09564

Item Number: CM0011553

Introduction

-conjugated rabbit monoclonal antibody detecting human CD7 (GP40, TP41) for flow cytometry. CD7 is a costimulatory transmembrane glycoprotein expressed on T and NK cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD7 mAb
Remarks/Alias	GP40; TP41; Tp40; LEU-9
Species	Human
Gene ID (Human)	924
Gene ID	924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag