

Rabbit Anti-Human Cd7 Monoclonal Antibody For Fc - P09564

Item Number: CM0005222

Introduction

Rabbit monoclonal antibody conjugated to , specific for human CD7. Validated for flow cytometry (FC). Recognizes an epitope within amino acids 26-180. Ideal for T-cell and NK cell immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD7 mAb
Remarks/Alias	GP40; TP41; Tp40; LEU-9
Species	Human
Gene ID (Human)	924
Gene ID	924
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag