

Rabbit Anti-Human Cd7 Monoclonal Antibody For Flow Cytometry - P09564

Item Number: CM0008669



Introduction

-conjugated rabbit monoclonal antibody against human CD7 for flow cytometry. Recognizes T-cell antigen CD7 (GP40/LEU-9). Suitable for T-cell and NK cell immunophenotyping. UniProt P09564.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD7 mAb
Remarks/Alias	CD7; GP40; LEU-9; TP41; Tp40; CD7 molecule
Species	Human
Gene ID (Human)	924
Gene ID	924
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag