


Rabbit Anti-Human Cd70 Monoclonal Antibody For Flow Cytometry - P32970

Item Number: CM0011864



Introduction

rabbit anti-human CD70 monoclonal antibody validated for flow cytometry. Recognizes CD70 (CD27L/TNFSF7), a T cell activation ligand, with human reactivity. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD70 mAb
Remarks/Alias	CD27L; LPFS3; CD27-L; CD27LG; TNFSF7; TNLG8A
Species	Human
GeneID (Human)	970
GeneID	970
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32970
Protein Weight	21kDa/23kDa
Shipping	Ice bag