


Rabbit Anti-Human Cd71 Monoclonal Antibody For Flow Cytometry - P02786

Item Number: CM0009738



Introduction

-conjugated rabbit monoclonal antibody targeting human CD71 (transferrin receptor 1) for flow cytometry. Recognizes human CD71, ~85kDa. Ideal for studies on iron uptake, receptor trafficking, and viral entry mechanisms.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD71 mAb
Remarks/Alias	T9; TR; TFR; p90; CD71; TFR1; TRFR; IMD46
Species	Human
Gene ID (Human)	7037
Gene ID	7037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human CD71 (NP_003225.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02786
Protein Weight	85kDa
Shipping	Ice bag