

Rabbit Anti-Human Cd72 Monoclonal Antibody For Fc - P21854

Item Number: CM0012140

Introduction

-conjugated rabbit monoclonal antibody targeting human CD72 (Lyb-2), validated for flow cytometry. This directly labeled antibody eliminates the need for secondary reagents, streamlining multi-color panel design for lymphocyte subset analysis. Suitable for B-cell research and immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD72 mAb
Remarks/Alias	CD72; CD72b; LYB2; CD72 molecule
Species	Human
GeneID (Human)	971
GeneID	971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 117-359 of human CD72(NP_001773.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21854
Protein Weight	40kDa
Shipping	Ice bag