

Rabbit Anti-Human Cd72 Monoclonal Antibody For Flow Cytometry - P21854

Item Number: CM0005243

Introduction

-conjugated rabbit monoclonal antibody against human CD72 (B-cell differentiation antigen) validated for flow cytometry. Ideal for studying BCR signaling, regulatory B/T cell interactions, and B-cell development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD72 mAb
Remarks/Alias	CD72; CD72b; LYB2; CD72 molecule
Species	Human
GenelD (Human)	971
GenelD	971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 117-359 of human CD72 (NP_001773.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21854
Protein Weight	40kDa
Shipping	Ice bag