

Rabbit Anti-Human Cd74 Monoclonal Antibody (Pe) For Flow Cytometry - P04233-2

Item Number: CM0009756

Introduction

Rabbit monoclonal antibody against human CD74, PE-conjugated, validated for flow cytometry. Recognizes the invariant chain of MHC class II, with multiple isoform weights from 18 to 34 kDa. Ideal for immunophenotyping and antigen presentation studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD74 mAb
Remarks/Alias	II; p33; CLIP; DHLAG; HLADG; Ia-GAMMA
Species	Human
GeneID (Human)	972
GeneID	972
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-232 of human CD74 (NP_004346.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04233-2
Protein Weight	18kDa/24kDa/26kDa/31kDa/34kDa
Shipping	Ice bag