

Rabbit Anti-Human Cd75 Monoclonal Antibody For Fc - P15907

Item Number: CM0008518



Introduction

Highly specific -conjugated rabbit monoclonal antibody targeting human CD75/ST6GAL1. Validated for flow cytometry, ideal for lymphocyte immunophenotyping and B-cell surface antigen studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD75 mAb
Remarks/Alias	ST6N; CDw75; SIAT1; ST6GalI
Species	Human
GeneID (Human)	6480
GeneID	6480
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15907
Protein Weight	20kDa/47kDa
Shipping	Ice bag