

Rabbit Anti-Human Cd75 Monoclonal Antibody For Fc - P15907

Item Number: CM0008142



Introduction

-conjugated rabbit monoclonal antibody specific for human CD75 (ST6GAL1), validated for flow cytometry. Enables detection of B-cell surface antigen CDw75 and sialyltransferase activity in research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD75 mAb
Remarks/Alias	ST6N; CDw75; SIAT1; ST6GalI
Species	Human
Gene ID	6480 (Human)
Immunogen	Recombinant fusion protein containing amino acids 44-406 of human CD75 (NP_001340845.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15907
Protein Weight	20kDa/47kDa
Shipping	Ice bag