

Rabbit Anti-Human Cd79A Monoclonal Antibody For Flow Cytometry (Intracellular) - P11912

Item Number: CM0010153



Introduction

-conjugated rabbit monoclonal antibody targeting human CD79a (MB-1, Ig-alpha). Validated for intracellular flow cytometry. Reliable for B-cell receptor signaling and B-cell lineage studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79a mAb
Remarks/Alias	IGA; MB1; MB-1; IGAAlpha
Species	Human
Gene ID (Human)	973
Gene ID	973
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11912
Protein Weight	20kDa/25kDa
Shipping	Ice bag