

Rabbit Anti-Human Cd79B/Igβ Mab For Flow Cytometry - P40259

Item Number: CM0006105



Introduction

Rabbit anti-Human CD79B/Igβ monoclonal antibody for flow cytometry targeting CD79B in human B-cells. Essential for studying B-cell receptor signaling, antigen presentation, and immune response mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79B/Igβ mAb
Remarks/Alias	B29; IGB; AGM6; Igbeta
Species	Human
Gene ID (Human)	974
Gene ID	974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79b/Igβ (NP_000617.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40259
Protein Weight	26kDa
Shipping	Ice bag