

Rabbit Anti-Human Cd79B/Igβ Monoclonal Antibody () For Flow Cytometry - P40259

Item Number: CM0008437



Introduction

High-quality rabbit monoclonal antibody against human CD79B/Igβ, conjugated to for flow cytometry. Specifically detects the extracellular domain of CD79B, a key B-cell antigen receptor complex component. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79B/Igβ mAb
Remarks/Alias	B29; IGB; AGM6; Igbeta
Species	Human
GeneID (Human)	974
GeneID	974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79b/Igβ (NP_000617.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40259
Protein Weight	26kDa
Shipping	Ice bag