

Rabbit Anti-Human Cd79B/Igβ Monoclonal Antibody For If/Fc - P40259

Item Number: CM0011231



Introduction

-conjugated rabbit monoclonal antibody specific for human CD79b (Igβ), a B-cell antigen receptor complex subunit. Validated for immunofluorescence (IF/ICC) and flow cytometry (FC). Ideal for B-cell research and diagnostics development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79b/Igβ mAb
Remarks/Alias	CD79B; AGM6; B29; IGB; CD79b molecule
Species	Human
Gene ID (human)	974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79b/Igβ (NP_000617.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P40259
Protein Weight	14kDa/26kDa
Shipping	Ice bag