

Rabbit Anti-Human Cd79B/Igβ Monoclonal Antibody For If/Icc, Fc - P40259

Item Number: CM0010772



Introduction

-conjugated rabbit monoclonal antibody targeting human CD79b/Igβ (B29, AGM6). Validated for IF/ICC and FC. CD79b is a B-cell antigen receptor complex protein crucial for BCR signaling and B-cell development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD79b/Igβ mAb
Remarks/Alias	CD79B; AGM6; B29; IGB; CD79b molecule
Species	Human
GeneID	974
Immunogen	Recombinant protein of human CD79b/Igβ
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P40259
Protein Weight	14kDa/26kDa
Shipping	Ice bag