

## Rabbit Anti-Human Cd79B/Igβ Monoclonal Antibody For If/Icc, Fc - P40259

Item Number: CM0010497



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD79b/Igβ, validated for immunofluorescence and flow cytometry. Specifically detects B-cell receptor complex beta chain, suitable for B-cell signaling and immunology studies.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD79b/Igβ mAb                                                                                       |
| Remarks/Alias    | CD79B; AGM6; B29; IGB; CD79b molecule                                                                                 |
| Species          | Human                                                                                                                 |
| Gene ID          | 974                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79b/Igβ (NP_000617.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | IF/ICC, FC                                                                                                            |
| Cross Reactivity | Human                                                                                                                 |
| SWISS            | P40259                                                                                                                |
| Protein Weight   | 14kDa/26kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                               |