

Rabbit Anti-Human Cd80/B7-1 Mab For Flow Cytometry - P33681

Item Number: CM0007840



Introduction

Rabbit anti-human CD80/B7-1 monoclonal antibody, conjugated to , for flow cytometry. Targets the T-lymphocyte activation antigen CD80, a critical costimulatory molecule in immune cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD80/B7-1 mAb
Remarks/Alias	CD80; B7; B7-1; B7.1; BB1; CD28LG; CD28LG1; LAB7; CD80 molecule
Species	Human
Gene ID (Human)	941
Gene ID	941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-242 of human CD80/B7-1 (NP_005182.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33681
Protein Weight	18kDa/29kDa/33kDa
Shipping	Ice bag