

Rabbit Anti-Human Cd80/B7-1 Mab - P33681

Item Number: CM0004263



Introduction

Anti-human CD80/B7-1 rabbit monoclonal antibody (clone) for flow cytometry. Specific to human CD80 (SWISS: P33681), a key costimulatory molecule in T-cell activation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD80/B7-1 mAb
Remarks/Alias	B7; BB1; B7-1; B7.1; LAB7; CD28LG; CD28LG1
Species	Human
GeneID (human)	941
GeneID	941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-242 of human CD80/B7-1 (NP_005182.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33681
Protein Weight	18kDa/29kDa/33kDa
Shipping	Ice bag