

Rabbit Anti-Human Cd80/B7-1 Monoclonal Antibody For Flow Cytometry - P33681

Item Number: CM0006032



Introduction

-conjugated rabbit monoclonal antibody targeting human CD80 (B7-1), validated for flow cytometry. Detects a key costimulatory molecule on antigen-presenting cells; Swiss-Prot: P33681.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD80/B7-1 mAb
Remarks/Alias	CD80; B7; B7-1; B7.1; BB1; CD28LG; CD28LG1; LAB7; CD80 molecule
Species	Human
GeneID (Human)	941
GeneID	941
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33681
Protein Weight	18kDa/29kDa/33kDa
Shipping	Ice bag