

Rabbit Anti-Human Cd80/B7-1 Monoclonal Antibody For Flow Cytometry - P33681

Item Number: CM0005348



Introduction

High-quality rabbit anti-human CD80/B7-1 monoclonal antibody for flow cytometry. Targets extracellular domain, reactive with human. Ideal for T-cell activation, immune checkpoint, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD80/B7-1 mAb
Remarks/Alias	B7; BB1; B7-1; B7.1; LAB7; CD28LG; CD28LG1
Species	Human
GeneID (Human)	941
GeneID	941
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-242 of human CD80/B7-1 (NP_005182.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33681
Protein Weight	18kDa/29kDa/33kDa
Shipping	Ice bag