


Rabbit Anti-Human Cd81 Monoclonal Antibody For Fc - P60033

Item Number: CM0009627



Introduction

Rabbit anti-Human CD81 Monoclonal Antibody for Flow Cytometry. Targets human CD81, a tetraspanin involved in immune signaling and viral entry. Suitable for FC studies on human B cells, T cells, and hepatocytes.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD81 mAb
Remarks/Alias	S5.7; CVID6; TAPA1; TSPAN28
Species	Human
GeneID (Human)	975
GeneID	975
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 114-199 of human CD81 (NP_004347.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P60033
Protein Weight	26kDa
Shipping	Ice bag