

Rabbit Anti-Human Cd82 Monoclonal Antibody For Flow Cytometry - P27701

Item Number: CM0008373



Introduction

-conjugated Rabbit anti-Human CD82 monoclonal antibody validated for flow cytometry. Targets the metastasis suppressor KAI1 (TSPAN27), involved in cell adhesion, signaling, and immune modulation. Ideal for oncology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD82 mAb
Remarks/Alias	R2, 4F9, C33, IA4, ST6, GR15, KAI1, SAR2, TSPAN27
Species	Human
GeneID (Human)	3732
GeneID	3732
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-267 of human CD82 (NP_002222.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P27701
Protein Weight	27kDa/30kDa
Shipping	Ice bag