


Rabbit Anti-Human Cd84 Monoclonal Antibody For Fc - Q9Uib8

Item Number: CM0005757



Introduction

-conjugated rabbit monoclonal antibody targeting human CD84 (SLAMF5), validated for flow cytometry. Recombinant immunogen (aa 22-225) detects human CD84; suitable for immune cell studies including NK cell, T-cell, B-cell, and dendritic cell analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD84 mAb
Remarks/Alias	CD84; LY9B; SLAMF5; hCD84; mCD84; CD84 molecule
Species	Human
Gene ID (Human)	8832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-225 of human CD84 (NP_001171808.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UIB8
Protein Weight	24kDa/26kDa/30kDa/31kDa/36kDa/38kDa
Shipping	Ice bag