


Rabbit Anti-Human Cd84 Monoclonal Antibody For Fc - Q9Uib8

Item Number: CM0004257



Introduction

-conjugated rabbit monoclonal antibody targeting human CD84 (SLAMF5). Validated for flow cytometry (FC) applications in human samples. Ideal for immune cell profiling and CD84 function studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD84 mAb
Remarks/Alias	CD84; LY9B; SLAMF5; hCD84; mCD84; CD84 molecule
Species	Human
GeneID (Human)	8832
GeneID	8832
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-225 of human CD84 (NP_001171808.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UIB8
Protein Weight	24kDa/26kDa/30kDa/31kDa/36kDa/38kDa
Shipping	Ice bag