

Rabbit Anti-Human Cd85D/E/I/J Monoclonal Antibody For Flow Cytometry - Q8N423 / Q8N6C8 / O75019 / Q8Nhl6

Item Number: CM0010669



Introduction

-conjugated rabbit monoclonal antibody targeting human CD85d (LILRB2) and cross-reactive with CD85e/i/j. Validated for flow cytometry, this antibody is ideal for immune tolerance and MHC class I receptor research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD85d/e/i/j mAb
Remarks/Alias	ILT4, LIR2, CD85D, ILT-4, LIR-2, MIR10, MIR-10
Species	Human
GeneID	10288/11026/11024/10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-461 of human CD85d (NP_001074447.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry
Cross Reactivity	Human
SWISS	Q8N423/Q8N6C8/O75019/Q8NHL6
Protein Weight	65kDa
Shipping	Ice bag