

Rabbit Anti-Human Cd85D/Ilt4 Mab For Fc - Q8N423

Item Number: CM0009000



Introduction

-conjugated rabbit monoclonal antibody targeting human CD85d/ILT4, validated for flow cytometry. Recognizes a broad spectrum of MHC class I antigens. Ideal for immune tolerance, maternal-fetal interface, and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD85d/ILT4 mAb
Remarks/Alias	ILT4; LIR2; CD85D; ILT-4; LIR-2; MIR10; MIR-10
Species	Human
GeneID (Human)	10288
GeneID	10288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-461 of human CD85d (NP_001074447.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8N423
Protein Weight	65kDa
Shipping	Ice bag