

Rabbit Anti-Human Cd85D/Itt4 Monoclonal Antibody For Fc - Q8N423

Item Number: CM0004544



Introduction

-conjugated rabbit anti-human CD85d/ILT4 monoclonal antibody for flow cytometry. Recognizes human LILRB2, a receptor for class I MHC antigens involved in immune tolerance. Validated for FC, cross-reacts with human. 65kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD85d/ILT4 mAb
Remarks/Alias	ILT4; LIR2; CD85D; ILT-4; LIR-2; MIR10; MIR-10
Species	Human
GeneID (Human)	10288
GeneID	10288
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8N423
Protein Weight	65kDa
Shipping	Ice bag