

Rabbit Anti-Human Cd85I/J Monoclonal Antibody For Flow Cytometry - O75019 / Q8Nhl6

Item Number: CM0010127

Introduction

High-quality -conjugated rabbit monoclonal antibody against human CD85i/j (LILRA1/LILRB1), validated for flow cytometry. Targets immune receptors expressed on monocytes and B-cells, ideal for immunophenotyping and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD85i/j Monoclonal Antibody
Remarks/Alias	LIR6; CD85i; LIR-6; CD85j; ILT-2; ILT2; LIR-1; LIR1; MIR-7; MIR7; PIR-B; PIRB
Species	Human
GeneID (Human)	11024/10859
GeneID	11024/10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-461 of human CD85i (NP_006854.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75019/Q8NHL6
Protein Weight	53kDa
Shipping	Ice bag